

HP 3395/3396 Quick Reference Card

Chart Control

- START** Starts run
- [SHFT] SEQ START** Starts a sequence of runs
- PLOT** Starts plot without integration
- STOP** Stops run or plot
- CHT SP 1 ENTER** Sets chart speed to 1 cm/min
- ZERO 10 ENTER** Positions plot 10% away from left margin
- ATT 2↑ 5 ENTER** Sets attenuation to 5
- [CTRL] KA** Sets top of form for 66 line pg (USA)
- [CTRL] VA4** Sets top of form for 70 line pg
- [CTRL] LFF** Form feed
- [CTRL] A** Advances paper 1/8 line
- [CTRL] W** Prints using large font
- [CTRL] D** Prints using small font
- [CTRL] N** Presentation plot
- [CTRL] O** Draft plot
- TIME** Prints run time during a run
- ZERO ENTER** Zeros chart during a run
- [SHFT] OP() 1 ENTER** Specifies plot type and presentation or draft plot
- [SHFT] OP() 5 ENTER** Specifies form feeds before and after reports, perforation skips in plots and reports, and large or small font

Data Storage

A: = Optional floppy disk drive (HP 3396 only)
E: = Optional EPROM
H: = Optional Host memory
M: = Main memory

- DIR [space] M: ENTER** Lists files on M: disk
- A: ENTER** Changes default disk from M: disk to A: disk (Default disk is M: unless otherwise specified)
- FORMAT [space] M:, VOL1, 200 ENTER** Initializes internal disk in drive M: names it VOL1 and divides it into 200 files
- COPY [space] A:Q00332D7.BNC, M:DATA391 ENTER**
Copies data file A:Q00332D7.BNC to M:DATA391
- RENAME [space] M:DATA391.BNC, M:Q00332D7 ENTER**
Renames data file M:DATA391.BNC to M:Q00332D7
- PURGE [space] METHOD1 ENTER** Deletes file named METHOD1
- [SHFT] OP() 2 ENTER** Specify whether to store data; determine either bunched or raw signal data; and where to store

Timetable

- LIST TIME ENTER** Lists all current timetable entries
- TIME .5 CHT SP 7 ENTER** Sets chart speed to 7 cm/sec at 0.5 min into the run
- DEL TIME ENTER** Deletes entire timetable
- DEL TIME 1 ENTER** Deletes all entries for 1 min
- TIME 1 INTG() 8 ENTER** Turns on Start/Stop marks 1 min into run
- TIME 2 INTG() -8 ENTER** Turns off Start/Stop marks 2 min into run

[SHFT] = Press SHIFT key and hold down while pressing next key
[CTRL] = Press CTRL key and hold down while pressing next key
[space] = Press the spacebar and release it

Status/Calendar

LIST **LIST** Lists current run parameters

R E A **ENTER** Checks system readiness

S Y **ENTER** Lists current system configuration

LIST **[SHFT]** **OP()** **2** **ENTER**

Lists current status for Option 2

T I M E **8 : 3 5 : 3 0** **ENTER**

Sets the clock to 8:35:30 AM

D A T E **6 / 1 / 9 0** **ENTER**

Sets the calendar to June 1, 1990

Reports

[SHFT] **AREA%** Prints uncalibrated percent report for data in default file

[SHFT] **REPORT** Prints calibrated report

N O T E P A D **ENTER**

Adds notes to printed report

I D [sp] **LAB1** **ENTER**

Adds identifier to report header

I D " " **ENTER**

Deletes identifier from report header

[SHFT] **OP()** **4** **ENTER** Specifies local report suppression, report title, amount label, inclusion of uncalibrated peaks, and report format

[SHFT] **OP()** **5** **ENTER** Specifies postrun report storage, external printing, and listing of additional information after a report.

Methods

[SHFT] **PREP** **[SHFT]** **METH** Prepares a method step by step

[SHFT] **EDIT** **[SHFT]** **METH** Edits an existing method

LIST **[SHFT]** **METH** **M E T H 6 2 5** **ENTER** Lists METH625 from default disk

[SHFT] **STORE** **[SHFT]** **METH** **M : M E T H 6 2 5** **ENTER** Stores current method parameters under METH625 on disk M:

[SHFT] **LOAD** **[SHFT]** **METH** **M : M E T H 6 2 5** **ENTER** Loads METH625 into integrator from M: disk

DEL **[SHFT]** **METH** **M : M E T H 6 2 5** **ENTER** Deletes METH625 from M: disk

[SHFT] **OP()** **5** **ENTER** Specifies form feeds before and after reports, perforation skips in plots and reports, and large or small font

[SHFT] **OP()** **6** **ENTER** Sends commands to specified remote device

LIST **[SHFT]** **OP()** **ENTER** Lists all option dialogs

Sequences

Syntax same as for methods using **[SHFT]** **SEQ** key ; see above

[SHFT] **OP()** **7** **ENTER** Specifies use of sample table in manual run

External Events (HP 19405A) (HP 3396 only)

EXT() **1** **ENTER** Turns on external event 1

EXT() **-1** **ENTER** Turns off external event 1

TIME **1** **EXT()** **5** **ENTER** Turns external event 5 on at 1 min into the run

DEL **TIME** **1** **EXT()** **ENTER** Deletes external events programmed for 1 min into run

DEL **TIME** **EXT()** Deletes entire external events timetable

LIST **TIME** **EXT()** Lists external event timetable

TIME **1** **EXT()** **ENTER** Sets equilibration time to 1 min and activates auto reset function

[SHFT] = Press SHIF'T key and hold down while pressing next key

[CTRL] = Press CTRL key and hold down while pressing next key

[space] = Press the spacebar and release it

Calibration/Calculations

[SHFT] **PREP** [SHFT] **CALIB** **ENTER** Prepares calibration step by step

[SHFT] **EDIT** [SHFT] **CALIB** **ENTER** Edits existing calibration

LIST [SHFT] **CALIB** **C A L 6 2 5** **ENTER** Lists calibration file CAL625 from default disk

[SHFT] **STORE** [SHFT] **CALIB** **M : C A L 6 2 5** **ENTER**
Stores current calibration parameters to CAL625 on disk M

[SHFT] **LOAD** [SHFT] **CALIB** **M : C A L 6 2 5** **ENTER**
Loads CAL625.CAL into integrator from M: disk

DEL [SHFT] **CALIB** **C A L 6 2 5** **ENTER** Deletes CAL625.CAL

[SHFT] **OP()** **3** **ENTER** Specifies response factor for uncalibrated peaks, calibration fit, retention time updating, peak number, internal standard amount, sample amount, and multiplication factor

[SHFT] **OP()** **4** **ENTER** Specifies area or height for calculations

[SHFT] **OP()** **7** **ENTER** Specifies values for internal standard amount, sample amount, and multiplication factor without sample table

Recalibration

[SHFT] **CALIB** **2** **ENTER** Manually initiates averaging for response data for level 2 of current calibration

[SHFT] **CALIB** **-2** **ENTER** Replaces old data with new data for level 2 of current calibration

Integration/Reintegration

PK WD **. 0 1** **ENTER** Sets peak width to .01 min

THRSH **5** **ENTER** Sets threshold to 5

THRSH **ENTER** Measures noise and sets threshold

THRSH **-** **ENTER** Aborts noise measurement

AR REJ **8 0 0** **ENTER** Sets area rejection limit to 800

TIME **. 0 1** **INTG()** **0** **ENTER** Sets baseline at .01 min

- 1** Sets baseline at next valley
- 2** Sets baseline at all valleys
- 3** Processes next peak as solvent peak
- 4** Turns off automatic solvent detection
- 5** Draws horizontal line
- 6** Measures and updates threshold
- 7** Turns off retention time labeling
- 8** Turns on Start/Stop marks
- 9** Turns off integration
- 10** Increments threshold
- 11** Inverts negative peaks
- 12** Clamps negative peaks
- 13** Shows functions 11 and 12
- 14** Starts peak sum window

A N **ENTER** Reintegrates data in default file

A N , I **ENTER** Reintegrates data in default file using original PK WD profile

A N [space] **D A T A 1** **ENTER** Reintegrates data in DATA1 file

A N [space] **D A T A 1 , I** **ENTER** Reintegrates data in DATA1 file using original PK WD profile

[SHFT] **OP()** **1** **ENTER** Specifies integration plot type

[SHFT] = Press SHIFT key and hold down while pressing next key

[CTRL] = Press CTRL key and hold down while pressing next key

[space] = Press the spacebar and release it

Applications Programs

At asterisk prompt (*) press appropriate key to start program.

1 File Manager

Manage your disk files with wildcards using the **COPY**, **DELETE**, **DIRECTORY**, **FORMAT**, **LIST**, and **RENAME** commands.

2 Batch Reprocess

Reprocesses existing data files using new method or sequence parameters. Calibration files can be reprocessed to update their calibration information.

- Updates calibration information of the current method.
- Updates calibration information of a specified calib file.

3 Plot Calibration Curve

Plots the response curves of calibrated peaks from a method or multi-level calibration file.

- Plots with defaults or enter plotting parameters.
- Plots selected peaks or the entire file.

5 Bar Coded Methods (HP 3396 Dual Channel only)

Automates runs without a sequence. The method, inj volume, cal level, and the number of injections is coded on each bar coded vial.

- Prepare bar code labels ahead of time.

6 Sequence Chaining

Chains a set of sequences together.

- Each sequence can be assigned an optional Autoscheduler file for scheduling postrun programs.

8 Plot Baseline

Replots the original chromatogram with its baseline.

- Set OP() 2 to store signal data and processed peak files.
- Store bunched data when possible.
- E:BASELINE.BAS can be scheduled as a postrun program.

9 Autoscheduler

Schedules postrun programs and provides access to the **AUTONAME** and **AUTO_2CH** (dual channel) programs.

- Press 0 to start an autoscheduled run or sequence.

Autoname

E:AUTONAME.BAS is an application program that automatically renames the signal data, processed peak, and report files after each run.

- Enter the file name prefix in the Autoscheduler dialog.
Prefix = TEST and Run#=002, then TEST002.BNC
Prefix = TEST*, Vial#=20, and Inj#=02, then TEST2002.BNC
If no prefix, sample names are read from the active sample table
- Schedule **E:AUTONAME.BAS** as a postrun program.
- Press 0 to start an autoscheduled run or sequence.

HP 5890 Dual Channel (HP 3396 only)

E:AUTO_2CH.BAS is the Dual Channel application program.

- The Dual Channel program must be the first postrun program on the buffered channel.
- Press 0 to start an autoscheduled run or sequence.

0 Auto Start

Start key for autoscheduled run or sequence.

