


Installation/Service Manual Economy Oven

Models 14EG (6520/6521), **25EG** (6522/6523), **45EG** (6524/6525),
25EM (6526/6527), **45EM** (6528/6529)



Thermo Electron Corporation
Millcreek Road, Box 649
Marietta, Ohio 45750
USA
Phone: 740-373-4763
Toll Free: 800-848-3080
FAX: 740-373-4189

Manual P/N3177886
Rev E Dated 11AUG05

NOTICE

THE MATERIAL IN THIS MANUAL IS FOR INFORMATION PURPOSES ONLY. THE CONTENTS AND THE PRODUCT IT DESCRIBES ARE SUBJECT TO CHANGE WITHOUT NOTICE. THERMO ELECTRON CORPORATION MAKES NO REPRESENTATIONS OR WARRANTIES WITH RESPECT TO THIS MANUAL. IN NO EVENT SHALL THERMO BE LIABLE FOR ANY DAMAGES, DIRECT OR INCIDENTAL, ARISING OUT OF OR RELATED TO THE USE OF THIS MANUAL.

For repair information or replacement parts assistance from the manufacturer, call Technical Services using our toll free telephone number.

888-213-1790
(FAX) 740-373-4189

REVISION STATUS

INDEX	DATE	AMENDED PAGES	NOTES
A	6/99		Initial release
B	11/01	6	Add "plastics can melt" note
C	11/02	2,5,7,8,11	High limit cutout for models 25 & 45 replaced by hydraulic thermostat (ECO 020809A)
D	9/04	12, 13	Corrected safety thermostat callout
E	8/05	36100121	new manual #, manufacture location

Contents

1. INTRODUCTION	1
2. UNPACKING AND DAMAGE	1
3. GENERAL INFORMATION	2
4. SPECIFICATIONS	3
5. INSTALLATION	4
6. EXPLANATION OF CONTROLS	4
7. OPERATION	5
8. MAINTENANCE & SERVICING	7
9. PARTS REPLACEMENT	8
10. ELECTRICAL DRAWINGS	11
11. REPLACEMENT PARTS LIST	13
12. WARRANTY	14



Warnings and Notices

Important: READ THESE NOTICES BEFORE INSTALLING THE UNIT

General Notices

AS A ROUTINE LABORATORY PRECAUTION, ALWAYS WEAR SAFETY GLASSES WHEN WORKING WITH THIS APPARATUS.

INSTALLATION SHOULD BE COMPLETED BY QUALIFIED PERSONNEL ONLY.

WITHIN THE MECHANICAL CONVECTION CABINETS, OBJECTS SHOULD NOT BE PLACED ON THE SHELVES IN SUCH A MANNER TO BLOCK THE MOVEMENT OF HEATED AIR INTO THE CHAMBER.

Electrical Warnings

FOR PERSONAL SAFETY, THIS APPARATUS MUST BE PROPERLY GROUNDED.

BE SURE THAT THE POWER SUPPLY IS OF THE SAME VOLTAGE AS SPECIFIED ON THE NAMEPLATE.

DANGEROUS VOLTAGES EXIST WITHIN THIS UNIT. SERVICE SHOULD BE PERFORMED ONLY BY QUALIFIED PERSONNEL. DISCONNECT THE UNIT FROM ITS ELECTRICAL SOURCE. REMOVE THE SHELVES AND THERMOMETER IF INSTALLED. DISCONNECTING ANY COMPONENT FROM THE CIRCUIT WITHOUT PRIOR REMOVAL OF THE POWER SOURCE MAY CAUSE DAMAGE TO OTHER CIRCUIT COMPONENTS.

DO NOT CRIMP OR SHARPLY BEND CAPILLARIES. PROVIDE ADEQUATE CLEARANCE BETWEEN CAPILLARY AND HEATER COILS TO PREVENT ELECTRICAL SHORT.

DISCONNECT OVEN FROM ITS POWER SOURCE BEFORE SERVICING.

1. INTRODUCTION

1.01 Your satisfaction and safety are important to Thermo and a complete understanding of this unit is necessary to attain these objectives.

1.02 As the user of this apparatus, you have the responsibility to understand the proper function and operational characteristics of it. This instruction manual should be thoroughly read and all operators given adequate training before attempting to place this unit in service. Awareness of the stated cautions and warnings, and compliance with recommended operating parameters -- together with maintenance requirements -- are important for safe and satisfactory operation. The unit should be used for its intended application; alterations or modifications will **VOID THE WARRANTY.**

WARNING

AS A ROUTINE LABORATORY PRECAUTION, ALWAYS WEAR SAFETY GLASSES WHEN WORKING WITH THIS APPARATUS.

1.03 This product is not intended, nor can it be used, as a sterile or patient connected device. In addition, this apparatus is not designed for use in Class I, II or III locations as defined by the National Electrical Code.

2. UNPACKING AND DAMAGE

2.01 Save all packing material if apparatus is received damaged. This merchandise was carefully packed and thoroughly inspected before leaving our factory.

2.02 Responsibility for safe delivery was assumed by the carrier upon acceptance of the shipment; therefore, claims for loss or damage sustained in transit must be made upon the carrier by the recipient as follows:

1. **Visible Loss or Damage:** Note any external evidence of loss or damage on the freight bill, or express receipt, and have it signed by the carrier's agent. Failure to adequately describe such external evidence of loss or damage may result in the carrier's refusing to honor your damaged claim. The form required to file such claim will be supplied by the carrier.
2. **Concealed Loss or Damage:** Concealed loss or damage means loss or damage which does not become apparent until the merchandise has been unpacked and inspected. Should either occur, make a written request for inspection by carrier's agent within fifteen (15) days of the delivery date; then file a claim with the carrier since the damage is the carrier's responsibility.

2.03 If you follow the above instructions carefully, we will guarantee our full support of your claim to be compensated for loss or concealed damage.

DO NOT -- FOR ANY REASON -- RETURN THIS UNIT WITHOUT FIRST OBTAINING AUTHORIZATION. In any correspondence to **Thermo**, please supply the nameplate data, including catalog number and serial number.

3. GENERAL INFORMATION

3.01 This instruction manual encompasses the following models and their specific electrical characteristics.

3.02 Models 14EG, 25EG and 45EG employ gravity convection as a method of heat transfer. Gravity convection is defined as the natural tendency for heated air to rise due to its change in density and mass.

3.02.1 Air is drawn into the chamber through openings in the bottom of the oven and then heated as it passes over the electric heating coils and up through the openings in the diffuser panel located on the bottom of the inner chamber. A limited amount of heated air is exhausted out of the top of the chamber through the vent shutter cap and the remaining air recirculates within the chamber.

3.02.2 Gravity convection ovens are ideal where forced air circulation cannot be tolerated and for situations demanding gentle curing or long term sample storage under closely controlled conditions. Gravity convection ovens are ideal for drying powders, soil samples, paper goods, semiconductors and cosmetics.

3.03 Models 25EM and 45EM use mechanical convection as a method of heat transfer. Mechanical convection can be defined as a positive and planned directional air flow or forced air circulation within the chamber.

3.03.1 Air is drawn into the chamber through vent tubes located in the bottom of the oven and is heated as it passes over the electric heating coils. The air is blown through openings in the bottom of the inner chamber. A limited amount of heated air is exhausted out the top of the chamber through the vent shutter cap and the remaining air recirculates within the chamber.

3.03.2 Mechanical convection ovens provide the most efficient means of heat transfer as well as the most reproducible test conditions for repeat operations. Mechanical convection allows for rapid heat up time for high density loads, shortened recovery periods after door openings and improved uniformity for extremely heat sensitive samples.

3.04 All models use a hydraulic thermostat for temperature safety limit control.

Economy Ovens

4. Specifications

Model Number	14EG	25EG	25EM	45EG	45EM
Convection Technique	Gravity	Gravity	Mechanical	Gravity	Mechanical
Temperature Control	Hydraulic	Electronic	Electronic	Electronic	Electronic
Temperature Display	0° to 250°C Thermometer	0° to 250°C Thermometer	0° to 250°C Thermometer	0° to 250°C Thermometer	0° to 250°C Thermometer
Temperature Range Sensitivity	65° to 210°C ±0.5°C	65° to 210°C ±0.3°C	50° to 210°C ±0.3°C	65° to 210°C ±0.3°C	50° to 210°C ±0.3°C
Dimensions (LxWxH) Chamber	14 x 12.5 x 13.75 in. 356 x 318 x 349 mm	15.5 x 18.5 x 15 in. 390 x 470 x 380 mm	15.5 x 18.5 x 15 in. 390 x 470 x 380 mm	15.5 x 18.5 x 27 in. 390 x 470 x 685 mm	15.5 x 18.5 x 27 in. 390 x 470 x 685 mm
Shelf (each)	13.5 x 12.5 in. 343 x 318 mm	15.5 x 18 in. 390 x 455 mm	15.5 x 18 in. 390 x 455 mm	15.5 x 18 in. 390 x 455 mm	15.5 x 18 in. 390 x 455 mm
Exterior*	19.5 x 18 x 28 in. 495 x 457 x 711 mm	21.5 x 24 x 28 in. 550 x 610 x 715 mm	21.5 x 24 x 28 in. 550 x 610 x 715 mm	21.5 x 24 x 40 in. 550 x 610 x 1020 mm	21.5 x 24 x 40 in. 550 x 610 x 1020 mm
Net Weight	69 lbs., 31 kg	93 lbs., 42.4 kg	103 lbs., 46.7 kg	120 lbs., 54.4 kg	130 lbs., 59 kg
Chamber Volume	1.4 cu.ft. 39 liters	2.5 cu.ft. 71.5 liters	2.5 cu.ft. 71.5 liters	4.5 cu.ft. 129 liters	4.5 cu.ft. 129 liters
Shelves** Supplied Maximum	2 6	2 6	2 6	2 10	2 10
Electrical Data*** (all models 50/60 Hz) 115V 230V	0.8 kW, 7.1 A 0.8 kW, 3.5 A	1.3 kW, 11.3 A 1.3 kW, 5.6 A	1.3 kW, 11.3 A 1.3 kW, 5.6 A	1.8 kW, 15.6 A 1.8 kW, 7.8 A	1.8 kW, 15.6 A 1.8 kW, 7.8 A
Shipping Data Dimensions (LxWxH)	24 x 23 x 30 in. 609 x 584 x 762 mm	25 x 26 x 31.5 in. 635 x 660 x 800 mm	25 x 26 x 31.5 in. 635 x 660 x 800 mm	25 x 26 x 43.5 in. 635 x 660 x 1105 mm	25 x 26 x 43.5 in. 635 x 660 x 1105 mm
Weight	79 lbs., 35.5kg	117 lbs., 53.1 kg	127 lbs., 59.6 kg	145 lbs., 65.8 kg	155 lbs., 70.3 kg
Volume	9.8 cu.ft. 0.28 cu. meters	11.7 cu.ft. 0.33 cu. meters	11.7 cu.ft. 0.33 cu. meters	16.3 cu.ft. 0.46 cu. meters	16.3 cu.ft. 0.46 cu. meters
Catalog Number 115 V 230 V	3166767 3166768	3166769 3166770	3166773 3166774	3166771 3166772	3166775 3166776

* Exterior height includes vent cap and adjustable feet

** Spacing between shelves is 2 inches (50mm)

*** All ovens equipped with line cord and plug

5. INSTALLATION

WARNING

INSTALLATION SHOULD BE COMPLETED BY QUALIFIED PERSONNEL ONLY.

5.01 Location - The most uniform operating conditions will be obtained by placing the oven in an area remote from drafts, ventilating outlets, radiators, and other rapidly changing ambient conditions. To assure proper ventilation allow a minimum of three (3) inches of clearance between the rear, top and sides of the oven and adjacent walls. If two or more ovens are to be placed side by side, then allow six (6) inches between them. The four legs on the bottom of the oven can be turned to raise or lower the corners so that it sits level on the table.

5.02 Electrical Connections -

1. The 14EG, 25EG and 25EM are provided with a 120V, 15 amp, three-prong (grounding) plug. The 45EG and 45EM come equipped with a 120V, 20 amp, three-prong (grounding) plug. Both plug configurations are industry standard, and mate with standard three-prong grounding wall receptacles to minimize the possibility of electric shock hazard from this apparatus. If in doubt the user should have the wall receptacle and circuit checked by a qualified electrician to make sure the receptacle can provide adequate current and is properly grounded.

WARNING

FOR PERSONAL SAFETY THIS APPARATUS MUST BE PROPERLY GROUNDED.

2. Where a standard two-prong wall receptacle is encountered, it is the personal responsibility and obligation of the user to have it replaced with a properly grounded three-prong wall receptacle. *DO NOT, UNDER ANY CIRCUMSTANCES, CUT OR REMOVE THE THIRD (GROUND) PRONG FROM THE POWER CORD. DO NOT USE A TWO-PRONG ADAPTER PLUG.*

5.03 Determine the total amount of current being used by other apparatus connected to the circuit that will be used for this apparatus. It is critical that the added current demand (see nameplate) of this and other equipment used on the same circuit does not exceed the rating of the fuse or circuit breaker.

CAUTION

BE SURE THAT THE POWER SUPPLY IS OF THE SAME VOLTAGE AS SPECIFIED ON THE NAMEPLATE.

6. EXPLANATION OF CONTROLS

6.01 On/Off Switch - The on/off switch controls the flow of all electric power to the oven. The blower motor in the mechanically convected models will always be in operation with the power switch ON.

6.02 Three-Heat Switch - Model 14EG has a three heat switch. This switch also functions as the main "on/off" switch. With this switch it is possible to supply maximum wattage to the unit for initial heat-up and then reduce the wattage when operating temperature is reached. Low setting may be used for temperatures up to approximately 115°C, medium setting from 115°C to 175°C and high setting above 175°C. Most uniform control is achieved when setting is just high enough to maintain desired temperature.

6.03 Temperature Control Knob - The Temperature Control knob is used to set the operating temperature. Numerical graduations do not refer to any specific temperature but are simply for reference. As you become familiar with operating your oven, record dial settings which correspond to your preferred operating temperatures.

6.04 Temperature Control Pilot Lamp - When the lamp above the Temperature Control is illuminated, this indicates electric current is being applied to the heater. It is normal for this lamp to cycle on and off during the operation of the oven.

6.05 Safety Control Knob - Model 14EG has a safety control knob to set the desired safety control temperature.

6.06 Safety Cutout Pilot Lamp - When the Safety Cutout Pilot Lamp is illuminated, this indicates that either the Temperature Control Thermostat or sensor has failed and that the oven is controlling from the Safety Controller.

6.07 Glass Thermometer - Used to indicate oven temperature and determine proper thermostat settings.




Figure 6.1 - Controls

7. OPERATION

7.01 Fully open the exhaust vent shutter cap on the top of the oven and keep it open at all times. However, if running the oven at the maximum rated temperature it may be necessary to turn the cap to a more closed position to retain heat.

Insert shelf supports into the holes punched in the side walls of the inner chamber. Insert the shelves into the shelf supports and try to keep an equal distance between shelves whenever possible. Never cover the shelves with foil or reduce their open surface area by more than 75%, this will greatly reduce convection and hence uniformity and control will suffer.

A thermometer is provided and should be inserted through the hole in the vent shutter cap such that the three metal fingers punched in the cap converge toward the center of the cap and the glass ring on the thermometer rests on the top of the cap. The metal fingers may have to be bent slightly to grip the thermometer.

1. Press the ON/OFF switch for power to the ON position. For the 14EG only, select the desired temperature range.
2. Rotate the safety control adjustment fully clockwise.
3. Rotate the Temperature Control knob to some median position.

4. Allow the chamber sufficient time to heat up and observe the temperature on the thermometer in the vent cap. The lamp above the Temperature Control knob should turn on and off at regular intervals which indicates the control has stabilized at a particular set point.

5. If the thermometer indicates that the chamber temperature is below the desired set point and the lamp above the Temperature Control has been turning on and off, rotate the Temperature Control slightly in a clockwise direction. Turn the Temperature Control slightly counterclockwise if the thermometer indicates the chamber temperature is above the desired set point.

6. Set the safety adjustment: When the thermometer indicates the chamber is at the desired temperature and the lamp above the temperature control is cycling on and off at fairly uniform intervals begin to turn the safety adjustment slowly counterclockwise. At the moment the safety lamp is illuminated, the lamp above the temperature control will turn off. Rotate clockwise until the lamp above the safety adjustment is off and the lamp above the temperature control is back on. The safety limit has now been set.

7.02 When you are finished using the oven, simply press the ON/OFF power switch to the OFF position. You may leave the temperature controls in their set positions if you wish to resume operation at that temperature some time later.

7.03 Loading

Although the gravity and mechanical convection ovens rely on different methods of air circulation, general loading procedures are applicable to both types and must be followed. It is important for uniformity and recovery that air circulation within the chamber is not restricted.

1. At least one (1) inch (2.5 cm) should be left between objects placed on the shelves.

NOTE: WITHIN THE MECHANICAL CONVECTION CABINETS, OBJECTS SHOULD NOT BE PLACED ON THE SHELVES IN SUCH A MANNER TO BLOCK THE MOVEMENT OF HEATED AIR INTO THE CHAMBER.

2. The bottom floor of the chamber must be kept free and clear of objects and never used as a shelf.

3. At no time should solid shelves be substituted for the shelves that are provided. Additional shelves and shelf supports are available from Precision.

7.04 Safety

1. DO NOT place any explosive, combustible or flammable materials in this chamber.

2. DO NOT place sealed containers in the chamber. Sealed containers, filled with materials, do not provide room from expansion or evaporation and can develop dangerous vapor pressure as the temperature increases.

3. Avoid placing plastic materials in the oven. Extreme temperatures may cause plastics to melt, posing a fire hazard. If plastics are placed in the oven, do not leave the oven unattended.

4. Avoid spillage of liquids or powders within the chamber. Clean all spills as soon as possible. Use caution if oven is still hot.

5. DO NOT evaporate noxious or poisonous fumes.

6. These ovens are not intended for food service or the preparation of meals.

7. DO NOT store containers filled with acidic or caustic solutions within the chambers, as vapors from these materials will attack the chamber interior and electrical components, thus voiding the warranty.

8. MAINTENANCE AND SERVICING

8.01 Troubleshooting

WARNING

DANGEROUS VOLTAGES EXIST WITHIN THIS UNIT. SERVICE SHOULD BE PERFORMED ONLY BY QUALIFIED PERSONNEL. DISCONNECT THE UNIT FROM ITS ELECTRICAL SOURCE. REMOVE THE SHELVES AND THERMOMETER IF INSTALLED. DISCONNECTING ANY COMPONENT FROM THE CIRCUIT WITHOUT PRIOR REMOVAL OF THE POWER SOURCE MAY CAUSE DAMAGE TO OTHER CIRCUIT COMPONENTS.

8.02 Temperature Variance or Fluctuation

1. Make sure vent shutter cap is not closed. Open to the maximum position.
2. Test unit when empty; if results are satisfactory the chamber was improperly loaded. Redistribute the load.
3. Be sure to allow ample time for an empty chamber to stabilize at the desired temperature setting. It could take over one hour to equilibrate depending upon the differences between ambient and operating temperatures. The mass of the load can also affect stabilization time.
4. Make certain severe line voltage fluctuations are not occurring.
5. Make certain all wire connections are secure at their terminals.
6. Make certain that an intermittent failure of the switch, thermostat or wiring has not occurred. Isolate the cause; repair or replace.

8.03 Heat Loss

Inspect door gasket to make certain it fits firmly against cabinet at all points. Replace gasket if torn or damaged.

Improper Door Closure - Inspect door latches which are spring loaded to see if they pull the door in tightly against the body of the oven. If the "finger" of the latch (in the door) has been sprung into the body of the latch, use a pen or similar object to pop it back out. Check the oven cabinet to see that it is level. Use a spirit level and turn the four adjustable feet on the bottom of the oven to make it level.

No Heat - If the chamber does not heat, first check the line voltage, circuit breakers and/or fuses of the line circuit. Check if all electrical connections are secure.

8.04 Heater Resistance Check

WARNING

DISCONNECT OVEN FROM ITS POWER SOURCE BEFORE PROCEEDING. REFER TO THE APPROPRIATE WIRING DIAGRAM AT THE END OF THE MANUAL AND LOCATE THE ELECTRICAL LEADS FOR THE HEATER. USE THE VALUES IN THE HEATER COLD RESISTANCE TABLE TO FIND THE APPROPRIATE VALUE. BE SURE TO DISCONNECT AT LEAST ONE HEATER LEAD FROM THE TERMINAL STRIP BEFORE TAKING THE RATING WITH YOUR OHM METER. AGAIN, IT IS IMPERATIVE THAT THE UNIT BE COMPLETELY DISCONNECTED FROM ITS ELECTRICAL POWER SOURCE BEFORE ANY READINGS ARE TAKEN.

If the heater is open (infinite resistance) it should be replaced. If the heater reads less than five ohms it is shorted and should be replaced. Check the resistance between each lead of the heater and a base metal point on the oven chassis. If there is less than infinite resistance (a million ohms or greater) between the heater and chassis, the heater is shorted to the ground and should be replaced. Also, inspect all wires leading to the heater for signs of shorting or electrical contact to chassis of oven.

9. PARTS REPLACEMENT

9.01 Heater Replacement

1. Disconnect power.
2. Remove the clamps securing the thermostat bulbs to the floor. Carefully bend the bulbs to a vertical position.
3. Remove two screws securing control panel and carefully lift control panel and lay it on it's face.
Model 14EG: Then remove three screws securing the front z-bracket, screws are located directly above the control panel. Remove one screw on back wall that secures the oven floor. Lift bottom floor straight up and out to remove it.
Model 25 & 45: Remove ten screws that secure floor and lift out.
4. Disconnect heater terminals and heater mounting screws. Carefully lift straight up and out.
5. Replace with new heater, reinstall in reverse order.

9.02 Probe Replacement, Models 25 & 45

1. Disconnect power.
2. Remove back panel screws.
3. Pull probe straight out.
4. Open control panel by removing top screws which secure it and disconnect probe terminals.
5. Replace with new probe, reinstall in reverse order.

9.03 Motor Replacement (Mechanical Models only)

1. Disconnect power.
2. Open control panel by removing top screws which secure it and disconnect motor and heater terminals.
3. Remove screws which secure lower floor and lift out.
4. Disconnect/remove heater mounting screws, remove heater and wires.
5. Detach blower wheel from motor shaft.

6. Turn oven on its back or side and remove bottom plate. Remove motor mounting bracket from underside of oven.
7. Replace with new motor, reinstall in reverse order.

9.04 Control Replacement, Models 25 & 45

1. Disconnect power
2. Open control panel by removing top screws which secure it and disconnect leads to control pcb.
3. Remove control pcb.
4. Replace with new control pcb, reinstall in reverse order.

9.05 Thermostat Replacement:

1. UNPLUG UNIT FROM ELECTRICAL SOURCE
2. Loosen thermostat bulbs by removing two clamps secured to chamber bottom floor.
3. Remove two screws securing control panel and carefully lift control panel and lay it on it's face.
4. Model 14 EG: Remove three screws securing the front z-bracket, screws are located directly above the control panel. Remove screw securing the oven floor. Lift bottom floor straight up and out, be careful with the thermostat bulbs.
Model 25 & 45: Remove ten screws that secure floor and lift out. Be careful with the thermostat bulb.
5. Remove knobs by loosening two set screws on each knob. (14 EG only)
6. Tag lead wires and detach from thermostats.

7. Model 14 EG: Remove screws securing thermostat to control panel. Remove two screws securing angle cover, located inside of inner chamber on the front left side. This will give you room to pull thermostat bulbs out.
Model 25 & 45: Remove thermostat bracket from control panel. Remove thermostat from bracket. Remove knob from thermostat.
8. Pull thermostat bulbs through chamber and out the front.
9. Retrace preceding steps to install the new thermostat assembly.

CAUTION

DO NOT CRIMP OR SHARPLY BEND CAPILLARIES. PROVIDE ADEQUATE CLEARANCE BETWEEN CAPILLARY AND HEATER COILS TO PREVENT ELECTRICAL SHORT.

9.06 Door and Latch Repair

A) Replace Door Assembly

1. Disconnect power
2. Remove door by unscrewing top hinge plate, be careful to hold door while removing hinge.
3. Lift door off lower hinge.
4. Replace door in reverse order.

B) Replace Door Latches


1. Replace door gasket by removing the screws which secure inner door liner. Reassemble in reverse order.
2. Replace door latches by removing inner door liner, door gasket and door insulation.
3. Unscrew door latch clips and replace latches.
4. Reassemble in reverse order.

C) Replace Door Handle

1. Replace door handle by removing inner door liner, door gasket and door insulation.
2. Remove screws which secure door handle and replace handle.
3. Reassemble in reverse order.

**9.07 Door Latch Grabber
(on body of unit)**


1. Disconnect power.
2. Remove screws which secure door latch grabber and remove them from front.
NOTE: Used retaining nuts will fall into oven body.
3. Insert replacement latch assembly through open slot and fasten with screws and washers provided.



INTERNAL SWITCH CONNECTIONS

OFF	N	→ NO CONNECTION
	LI	→ NO CONNECTION
LOW	N	→ 3
	LI	→ 2
MED	N	→ 3
	LI	→ 1
HIGH	N	→ 2 & 3
	LI	→ 1

WIRING DIAGRAM - MODEL 14 EG



HEATER COLD RESISTANCE BETWEEN TB1-4 AND TB1-2		
HEATER PART NO.	WATTS	OHMS
34247439	1100	47
34247440	1600	32.5

WIRING DIAGRAM
 MODELS 25EG, 25EM, 45EG, 45EM

REPLACEMENT PARTS LIST - ECONOMY OVEN - ALL MODELS			
	MODEL 14EG	MODEL 25 (ALL)	MODEL 45 (ALL)
VENT SHUTTER CAP	3177386		
DOOR HANDLE	3174923		
KNOB, TEMPERATURE CONTROL	3174903		
THERMOMETER 0-250°C	3175996		
TOP HINGE	316231		
BOTTOM HINGE	316232		
DOOR HINGE BUSHING	3177622		
SHELF KIT	3166212	3166214	
HEATER ASSY, 115V	3175450	3175511	3175512
HEATER ASSY, 220V	3175473	3175511	3175512
LATCH REPLACEMENT KIT	3167037		
SHELF SUPPORT	316239	316225	
WIRE SHELF	3177372	3177585	
DOOR GASKET	3174681		
SAFETY THERMOSTAT	3175234		
SWITCH, DPST, NON-LIGHTED	N/A	3175318	
SWITCH, 3 HEAT	3175269	N/A	
LEVELING FOOT, 5/16 - 18	3175336		
TEMPERATURE CONTROLLER	3175234	3176730	
TEMPERATURE SENSOR	N/A	3176735	
ECONOMY OVEN HARNESS	N/A	3176727	
TERM BLOCK 6 POSITION	N/A	3175416	
DOOR SUPPORT ANGLE	316236		
DOOR STRIKER ANGLE	316237		
LIGHT ASSEMBLY, 115V	3162435	3162433	
LIGHT ASSEMBLY, 230V	3162436	3162434	
KIT, OVEN PACKAGING	3167299	3167252	3167253
CORD SET 115V	3165481	3176552	3176562
CORD SET 220V	3165488	3176559	
MECHANICAL MODELS ONLY			
BLOWER WHEEL	N/A	3175897	
MOTOR ASSY, 115V	N/A	3166987	
MOTOR ASSY, 230V	N/A	3167015	

THERMO ELECTRON CORPORATION STANDARD PRODUCT WARRANTY

The Warranty Period starts two weeks from the date your equipment is shipped from our facility. This allows for shipping time so the warranty will go into effect at approximately the same time your equipment is delivered. The warranty protection extends to any subsequent owner during the first year warranty period.

During the first year, component parts proven to be non-conforming in materials or workmanship will be repaired or replaced at Thermo's expense, labor included. Installation and calibration are not covered by this warranty agreement. The Technical Services Department must be contacted for warranty determination and direction prior to performance of any repairs. Expendable items, glass, filters and gaskets are excluded from this warranty.

Replacement or repair of components parts or equipment under this warranty shall not extend the warranty to either the equipment or to the component part beyond the original warranty period. The Technical Services Department must give prior approval for return of any components or equipment. At Thermo's option, all non-conforming parts must be returned to Thermo Electron Corporation postage paid and replacement parts are shipped FOB destination.

THIS WARRANTY IS EXCLUSIVE AND IN LIEU OF ALL OTHER WARRANTIES, WHETHER WRITTEN, ORAL OR IMPLIED. NO WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE SHALL APPLY.

Thermo shall not be liable for any indirect or consequential damages including, without limitation, damages relating to lost profits or loss of products.

Your local Thermo Sales Office is ready to help with comprehensive site preparation information before your equipment arrives. Printed instruction manuals carefully detail equipment installation, operation and preventive maintenance.

If equipment service is required, please call your Technical Services Office at 1-888-213-1790 (USA and Canada) or 1-740-373-4763. We're ready to answer your questions on equipment warranty, operation, maintenance, service and special application. Outside the USA, contact your local distributor for warranty information.



